

HI Titers of the Ferret Sera

date of acceptance: 2019.7.15

item quantity 18

	NO.	Testing Date	ST/92/2009	A/California/07/2009	CK/HK/YU22/2002	ST/55/2012	B/Brisbane/60/2008	B/HK/461141/2010	ST/1163/2009	CK/HK/N T155/2008	G1	ST/411/2014	非特异反应	备注
			H1N1	H1N1	H5N1	H3N2	FluB	FluB	B2506(Malaysia)	H9N2	H9N2	H7N9		
1	8593	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
2	8594	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
3	8595	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
4	8596	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
5	8597	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
6	8598	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
7	8599	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
8	8600	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
9	8601	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
10	8602	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
11	8603	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
12	8604	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
13	8605	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
14	8606	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
15	8607	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
16	8608	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
17	8609	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
18	8610	2019. 7. 15	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)		
remarks			1. All sera were treated with RDE (receptor-destroying enzyme) and absorbed with Turkey Red Blood Cells to											
			2. The antibody titer was determined by testing serial two-fold dilutions from 1/20 to 1/2560 with haemagglutination											
			3. Serum with HI titer ≥ 40 can be regarded as positive for influenza virus.											

presented International Infection And Immunity Research Center, Shantou University

date: 2019-7-15

兽医VETERINARIAN: 侯学超

签名SIGNATURE:

签字日期 DATE OF ISSUE:

