

TD.84224 EFA Deficient Diet

Formula

Formula	g/Kg
Casein, "Vitamin-Free" Test	192.0
DL-Methionine	3.0
Sucrose	656.2
Coconut Oil, hydrogenated	50.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	3.8
Vitamin Mix, Teklad (40060)	10.0

Footnote

A diet with approximately 5% fat coming from hydrogenated coconut oil, which contains virtually no essential fatty acids. "Vitamin-Free" Test Casein (alcohol-extracted) is used to eliminate the small amount of fat that would come from regular casein.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	17.6	18.5
Carbohydrate	66.5	69.7
Fat	5.0	11.8

Kcal/g 3.8

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- Essential Fatty Acid (EFA)
- Hydrogenated Coconut Oil
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better