

# TD.92163 Selenium Deficient Diet

Formula	g/Kg
Torula Yeast	300.0
DL-Methionine	3.0
Sucrose	591.0
Corn Oil	50.0
Mineral Mix, Se Deficient (80313)	35.0
Calcium Carbonate	11.0
Vitamin Mix, Teklad (40060)	10.0

## Key Features

- Purified Diet
- Selenium
- Torula Yeast
- Rodent

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Footnote

Estimated background selenium for this kind of formula is less than 0.01 ppm Se. This is slightly lower than for casein type diets. TD.96363 is a possible ingredient matched control with selenium added at 0.1 mg/kg diet from sodium selenite. For questions about this or other selenium adjusted diets [askanutritionist@harlan.com](mailto:askanutritionist@harlan.com).

## Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	15.3	16.5
Carbohydrate	64.1	69.2
Fat	5.9	14.3

**Kcal/g 3.7**

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

*Teklad Diets are designed & manufactured for research purposes only.*

## Speak With A Nutritionist

- (800) 483-5523 *option 4*
- [askanutritionist@harlan.com](mailto:askanutritionist@harlan.com)

Harlan Laboratories, Inc.  
PO Box 44220  
Madison, WI 53744-4220



[www.harlan.com](http://www.harlan.com)

## Contact Us

- Obtain Pricing · Check Order Status ·
- [tekladinfo@harlan.com](mailto:tekladinfo@harlan.com)
- (800) 483-5523 *option 3*

## International Inquiry (Outside USA or Canada)

- [askanutritionist@harlan.com](mailto:askanutritionist@harlan.com)

## Place Your Order (USA & Canada)

- Please Choose One ·
- <https://webapps.harlan.com/teklad-orders/>
- [tekladinfo@harlan.com](mailto:tekladinfo@harlan.com)
- (800) 483-5523 *option 3*
- (608) 277-2066 *facsimile*

*Helping you do research better*