

CA.170481 AIN-76A Purified Diet (Rats/Mice)

Formula

Formula	g/Kg
Casein	200.0
DL-Methionine	3.0
Sucrose	499.99
Corn Starch	150.0
Corn Oil	50.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Vitamin Mix, AIN-76A (40077)	10.0
Choline Bitartrate	2.0
Ethoxyquin, antioxidant	0.01

Footnote

References: Second Report of ad hoc Committee on Standards for Nutritional Studies (1980). J. Nutrition 110, 1726. Report of the American Institute on Nutrition ad hoc Committee on Standards for Nutritional Studies (1977). J. Nutrition 107, 1340-1348.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	17.7	18.8
Carbohydrate	64.9	68.8
Fat	5.2	12.4

Kcal/g 3.8

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- American Institute of Nutrition
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation not advised
 - Contact a nutritionist for recommendations

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better