

TD.90016 6% Protein Diet

Formula

Formula	g/Kg
Casein	69.0
DL-Methionine	0.9
Sucrose	571.8
Corn Starch	200.0
Corn Oil	53.9
Cellulose	57.82
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.01
Mineral Mix, Ca-P Deficient (79055)	13.37
Calcium Phosphate, dibasic	21.6
Calcium Carbonate	1.6

Footnote

The 6% protein diet is in a popular series of adjusted protein diets. These diets are isocaloric and are matched for calcium (0.7%) and phosphorus (0.54%).

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	6.1	6.5
Carbohydrate	75.6	80.4
Fat	5.5	13.2

Kcal/g 3.8

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- Protein Adjusted Series
- Isocaloric
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

· Outside U.S.A. or Canada ·

- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

· Place Order · Obtain Pricing ·
· Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better