

CA.40060 Vitamin Mix, Teklad

Formula	g/Kg
p-Aminobenzoic Acid	11.0132
Vitamin C, ascorbic acid, coated (97.5%)	101.6604
Biotin	0.0441
Vitamin B ₁₂ (0.1% in mannitol)	2.9736
Calcium Pantothenate	6.6079
Choline Dihydrogen Citrate	349.6916
Folic Acid	0.1982
Inositol	11.0132
Vitamin K ₃ , menadione	4.9559
Niacin	9.9119
Pyridoxine HCl	2.2026
Riboflavin	2.2026
Thiamin HCl	2.2026
Vitamin A Palmitate (500,000 IU/g)	3.9648
Vitamin D ₃ , cholecalciferol (500,000 IU/g)	0.4405
Vitamin E, DL-alpha tocopheryl acetate (500 IU/g)	24.2291
Corn Starch	466.6878

Footnote

Designed as a general purpose, well-fortified vitamin mix. The recommended use level is 1% in rodent diets, but may be used in diets for other species as well. This mix includes a source of choline and vitamin C.

Teklad Mixes are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com



Harlan, Harlan Laboratories, Helping you do research better,
and the Harlan logo are trademarks and trade names of Harlan Laboratories, Inc.
© 2008 Harlan Laboratories, Inc.

Key Features

- Vitamin Mix
- General Purpose
- Commonly Used

Key Planning Information

- Products are made fresh to order
- Lead time 2 weeks
- Store product at 4°C or lower
- Use within 6 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number



Product Specific Information

- Powder
- Minimum order 500 gms

Options (Fees Will Apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·
- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better