

TD.94048 AIN-93M Purified Diet

Formula

| Formula | g/Kg |
|---------------------------------|---------|
| Casein | 140.0 |
| L-Cystine | 1.8 |
| Corn Starch | 465.692 |
| Maltodextrin | 155.0 |
| Sucrose | 100.0 |
| Soybean Oil | 40.0 |
| Cellulose | 50.0 |
| Mineral Mix, AIN-93M-MX (94049) | 35.0 |
| Vitamin Mix, AIN-93-VX (94047) | 10.0 |
| Choline Bitartrate | 2.5 |
| TBHQ, antioxidant | 0.008 |

Footnote

Reference: J. Nutr. 123:1939-1951, 1993. Formulated for the maintenance of adult, non-reproducing rodents.

Selected Nutrient Information¹

| | % by weight | % kcal from |
|--------------|-------------|-------------|
| Protein | 12.4 | 13.8 |
| Carbohydrate | 68.4 | 76.0 |
| Fat | 4.1 | 10.2 |

Kcal/g 3.6

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- American Institute of Nutrition
- Maintenance
- Commonly Used

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation not advised
 - Contact a nutritionist for recommendations

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better