

TD.86143 Vitamin A Deficient Diet

Formula

	g/Kg
Casein, "Vitamin-Free" Test	193.0
DL-Methionine	3.0
Sucrose	506.3776
Corn Starch	150.0
Cottonseed Oil	50.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, w/o choline, A, D, E (83171)	5.0
Choline Dihydrogen Citrate	3.497
Vitamin D3, cholecalciferol (500,000 IU/g)	0.0044
Vitamin E, DL-alpha tocopheryl acetate (1000 IU/g)	0.121

Footnote

"Vitamin-Free" Test Casein (alcohol-extracted) is used to further limit vitamin contribution from the protein source. TD.91280 is a possible ingredient matched control with vitamin A added back at 20 IU/g. For questions about this or other vitamin A adjusted diets askanutritionist@harlan.com.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	17.7	18.8
Carbohydrate	64.9	69.2
Fat	5.0	12.0

Kcal/g 3.8

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523 option 4
- askanutritionist@harlan.com

Harlan Laboratories, Inc.
PO Box 44220
Madison, WI 53744-4220



www.harlan.com

Key Features

- Purified Diet
- Vitamin A
- Alcohol-Extracted Casein
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Contact Us

- Obtain Pricing · Check Order Status ·
- tekladinfo@harlan.com
- (800) 483-5523 option 3

International Inquiry (Outside USA or Canada)

- askanutritionist@harlan.com

Place Your Order (USA & Canada)

- Please Choose One ·
- https://webapps.harlan.com/teklad-orders/
- tekladinfo@harlan.com
- (800) 483-5523 option 3
- (608) 277-2066 facsimile

Helping you do research better