TD.96208 0.49% NaCl Diet

Form	nula	g/Kg
Wheat		350.0
Corn		308.79
Soybea	an Meal (48%)	190.8
Corn G	luten Meal (60%)	50.0
Alfalfa	Meal (17%), dehydrated	30.0
Corn O	il	33.0
Dicalci	um Phosphate, FG (18.5% P, 21% Ca)	14.0
Calciur	n Carbonate, FG (38%)	12.0
Sodium	n Chloride	4.9
Minera	Mix, TSD (80318)	1.5
Vitamir	Mix, TSD (81125)	3.0
DL-Met	thionine, FG (99%)	1.0
L-Lysin	e HCl, FG (78%)	1.0
Ethoxy	quin, antioxidant	0.01

Key Features

- Natural Ingredient Diet
- Adjusted NaCl Series
- Rodent

Key Planning Information

- · Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Footnote

This diet is a modification of TD.90228. This formula is part of a series of adjusted NaCl formulas of similar composition. Approximate nutrient content: 19% protein, 5% fat, 0.9% Ca, 0.65% P, 0.15% Mg, 0.8% K, 0.2% Na, and 0.36% Cl. Contact Technical Services for details.

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Selected Nutrient Information ¹			
	% by weight	% kcal from	
Protein	19.3	23.8	
Carbohydrate ²	50.2	61.8	
Fat	5.2	14.4	

Kcal/g 3.2

¹Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

(800) 483-5523

www.harlan.com

askanutritionist@harlan.com



·Outside U.S.A or Canada ·

International Inquiry

- askanutritionist@harlan.com

Place Your Order (U.S.A & Canada)

- · Place Order · Obtain Pricing ·
 - · Check Order Status ·
- (800) 483-5523
- (608) 277-2066 facsimile
- tekladinfo@harlan.com

Helping you do research better

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220