

TD.88232 Basal Mix (without 20% Oil)

Formula

Formula	g/Kg
Casein	300.0
DL-Methionine	4.5
Sucrose	373.5
Corn Starch	187.5
Cellulose	62.5
Mineral Mix, AIN-76 (170915)	52.5
Calcium Carbonate	4.5
Vitamin Mix, Teklad (40060)	15.0

Footnote

This basal mix is designed to be used at the rate of 800 g/Kg of diet, in conjunction with 200 g of the oil source(s).

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	26.6	32.1
Carbohydrate	55.5	67.0
Fat	0.3	0.8

Kcal/g 3.3

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Basal Mix
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

· Outside U.S.A. or Canada ·

- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

· Place Order · Obtain Pricing ·
· Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better