

# TD.90018 40% Protein Diet

## Formula

Formula	g/Kg
Casein	460.0
Sucrose	231.82
Corn Starch	200.0
Corn Oil	50.0
Cellulose	15.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.01
Mineral Mix, Ca-P Deficient (79055)	13.37
Calcium Phosphate, dibasic	9.6
Calcium Carbonate	10.2

## Footnote

The 40% protein diet in a popular series of adjusted protein diets. These diets are isocaloric and are matched for calcium (0.7%) and phosphorus (0.54%).

## Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	40.0	42.6
Carbohydrate	41.6	44.3
Fat	5.5	13.2

**Kcal/g 3.8**

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

*Teklad Diets are designed & manufactured for research purposes only.*

## Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

[www.harlan.com](http://www.harlan.com)

## Key Features

- Purified Diet
- Protein Adjusted Series
- Isocaloric
- Rodent

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

## Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



*Helping you do research better*