

TD.06416 Adjusted Calories Diet (10/Fat)

Formula

	g/Kg
Casein	210.0
L-Cystine	3.0
Corn Starch	280.0
Maltodextrin	50.0
Sucrose	325.0
Lard	20.0
Soybean Oil	20.0
Cellulose	37.15
Mineral Mix, AIN-93G-MX (94046)	35.0
Calcium Phosphate, dibasic	2.0
Vitamin Mix, AIN-93-VX (94047)	15.0
Choline Bitartrate	2.75
Yellow Food Color	0.1

Footnote

Approx. 10% of total calories come from fat. Designed with similarities to Research Diets, Inc. formula D12450B. For the series TD.06414-TD.06416. Approximate fatty acid profile (% of total fat): 29% saturated, 37% monounsaturated, 34% polyunsaturated. Modifications with lower sucrose include TD.08806 with 9% and TD.110675 with 16%. Contact us at askanutritionist@harlan.com for other options.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	18.6	20.1
Carbohydrate	64.7	69.8
Fat	4.2	10.2

Kcal/g 3.7

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com



Key Features

- Purified Diet
- Diet Induced Obesity
- Control

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

· Outside U.S.A. or Canada ·

- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

· Place Order · Obtain Pricing ·

· Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com

Helping you do research better