

# TD.03584 35% Lard Diet (Adj., No C)

## Formula

| Formula                                   | g/Kg   |
|---|--------|
| Casein                                    | 230.0  |
| DL-Methionine                             | 3.5    |
| Sucrose                                   | 150.0  |
| Maltodextrin                              | 191.35 |
| Lard                                      | 350.0  |
| Mineral Mix, AIN-93G-MX (94046)           | 50.0   |
| Calcium Phosphate, monobasic, monohydrate | 2.5    |
| Vitamin Mix, AIN-93-VX (94047)            | 18.4   |
| Choline Bitartrate                        | 4.25   |

## Footnote

A diet used to induce obesity. A modification of TD.03307, using a vitamin mix without vitamin C. Approximate fatty acid profile (% total fat): 40% saturated, 50% monounsaturated, 10% polyunsaturated. Cited in Circ Res (2005) 96:1178-1184.

## Selected Nutrient Information<sup>1</sup>

|              | % by weight | % kcal from |
|--------------|-------------|-------------|
| Protein      | 20.4        | 15.0        |
| Carbohydrate | 36.1        | 26.6        |
| Fat          | 35.2        | 58.4        |

**Kcal/g 5.4**

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

*Teklad Diets are designed & manufactured for research purposes only.*

## Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

[www.harlan.com](http://www.harlan.com)

## Key Features

- Purified Diet
- Diet Induced Obesity
- Lard
- Dough Form or Pellets

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## International Inquiry

· Outside U.S.A. or Canada ·

- askanutritionist@harlan.com

## Place Your Order (U.S.A. & Canada)

· Place Order · Obtain Pricing ·  
· Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



*Helping you do research better*