

CA.170915 Mineral Mix AIN-76

Formula

	g/Kg
Sucrose, fine ground	118.03
Calcium Phosphate, dibasic	500.0
Sodium Chloride	74.0
Potassium Citrate, monohydrate	220.0
Potassium Sulfate	52.0
Magnesium Oxide	24.0
Manganous Carbonate	3.5
Ferric Citrate	6.0
Zinc Carbonate	1.6
Cupric Carbonate	0.3
Potassium Iodate	0.01
Sodium Selenite, pentahydrate	0.01
Chromium Potassium Sulfate, dodecahydrate	0.55

Footnote

Designed to be used at 3.5% of the diets. Reference: Report of the American Institute on Nutrition Ad hoc Committee on Standards for Nutritional Studies (1977). J. Nutrition 107, 1340-1348.

Teklad Mixes are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com



Key Features

- Mineral Mix
- AIN-76A
- Commonly Used
- Rodent

Key Planning Information

- Products are made fresh to order
- Lead time 2 weeks
- Store in cool, dry conditions
- Use within 12 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number



Product Specific Information

- Powder
- Minimum order 500 gms

Options (Fees Will Apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·
- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com

Helping you do research better