

TD.94045 AIN-93G Purified Diet

Formula

Formula	g/Kg
Casein	200.0
L-Cystine	3.0
Corn Starch	397.486
Maltodextrin	132.0
Sucrose	100.0
Soybean Oil	70.0
Cellulose	50.0
Mineral Mix, AIN-93G-MX (94046)	35.0
Vitamin Mix, AIN-93-VX (94047)	10.0
Choline Bitartrate	2.5
TBHQ, antioxidant	0.014

Footnote

Reference: J. Nutr. 123:1939-1951, 1993. Formulated for the growth, pregnancy, and lactational phases of rodents.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	17.7	18.8
Carbohydrate	60.1	63.9
Fat	7.2	17.2

Kcal/g 3.8

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- American Institute of Nutrition
- Growth
- Commonly Used

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation not advised
 - Contact a nutritionist for recommendations

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

· Outside U.S.A. or Canada ·

- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

· Place Order · Obtain Pricing ·
· Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better