

TD.98090 70% Carbohydrate Diet

Formula

Formula	g/Kg
Casein	200.0
DL-Methionine	3.0
Sucrose	645.6
Corn Starch	20.0
Maltodextrin	20.0
Soybean Oil	50.0
Cellulose	9.89
Vitamin Mix, AIN-93-VX (94047)	10.0
Choline Bitartrate	2.5
TBHQ, antioxidant	0.01
Mineral Mix, AIN-93G-MX (94046)	35.0
Calcium Phosphate, dibasic	4.0

Footnote

A diet with 70% carbohydrate as expressed by weight and % of total kcal. Sucrose makes up approximately 66% of the diet (includes contribution from premixes).

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	17.7	17.8
Carbohydrate	70.0	70.4
Fat	5.2	11.8

Kcal/g 4.0

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- Adjusted Carbohydrate
- Sucrose
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation not advised
 - Contact a nutritionist for recommendations

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better