

TD.00235 Basal Mix (93%)

Formula	g/Kg
Casein	215.06
L-Cystine	3.23
Corn Starch	427.3749
Maltodextrin	141.94
Sucrose	107.53
Cellulose	53.76
Mineral Mix, AIN-93G-MX (94046)	37.64
Vitamin Mix, AIN-93-VX (94047)	10.76
Choline Bitartrate	2.69
TBHQ, antioxidant	0.0151

Key Features

- Basal Mix
- AIN-93G Modification
- For Fat Addition
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- Powder
- Minimum order 3 Kg
- Irradiation not advised
 - Contact a nutritionist for recommendations

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·
- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better

Footnote

This basal mix is designed to be used at the rate of 930 g/kg diet plus fat at 70 g/kg to complete the diet. A modification from diet AIN-93G (TD.94045).

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	19.0	22.6
Carbohydrate	64.6	76.9
Fat	0.2	0.5

Kcal/g 3.4

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com