

# TD.01383 2018 + 2% Cholesterol

Formula	g/Kg
2018, Teklad Global 18% Protein Rodent Diet	980.0
Cholesterol	20.0

### Footnote

Adds 2% cholesterol to the 2018 Teklad Global 18% Protein Rodent Diet.

### Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	17.8	22.9
Carbohydrate <sup>2</sup>	47.0	60.5
Fat	5.7	16.5

**Kcal/g 3.1**

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

*Teklad Diets are designed & manufactured for research purposes only.*

### Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220  
[www.harlan.com](http://www.harlan.com)



## Key Features

- Standard Diet Base
- 2018
- Cholesterol
- Atherogenic

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### International Inquiry

- Outside U.S.A. or Canada ·
- askanutritionist@harlan.com

### Place Your Order (U.S.A. & Canada)

- Place Order · Obtain Pricing ·
- Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



*Helping you do research better*