

TD.93075 Adjusted Fat Diet (55/Fat)**Formula**

	g/Kg
Casein	289.0
DL-Methionine	3.33
Corn Starch	207.3
Sucrose	90.5
Vegetable Shortening, hydrogenated (Primex)	274.1
Corn Oil	16.0
Cellulose	53.12
Vitamin Mix, Teklad (40060)	13.33
Mineral Mix, AIN-76 (170915)	46.66
Calcium Carbonate	6.66

Footnote

A diet used to induce obesity. Approximate fatty acid profile (% total fat): 28% saturated, 30% trans, 28% monounsaturated (cis), 14% polyunsaturated (cis). Cited in Diabetes (2004) 53:3274-3285.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	25.5	21.2
Carbohydrate	28.8	24.0
Fat	29.3	54.8

Kcal/g 4.8

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220

www.harlan.com

Key Features

- Purified Diet
- Diet Induced Obesity
- Trans Fatty Acids
- Dough Form

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated

**Product Specific Information**

- Powder (dough)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)

International Inquiry

· Outside U.S.A. or Canada ·

- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

· Place Order · Obtain Pricing ·
· Check Order Status ·

- (800) 483-5523
- (608) 277-2066 *facsimile*
- tekladinfo@harlan.com



Helping you do research better